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Use and Maintenance Instructions for

Torque Ratchet (Art.No. 5012303002)

Document No. 5014040111 - Revision 06/2016

1 Using the Torqu Ratchet

The torque ratchet consists of the ratchet head (Fig. 1-b), the ratchet handle (Fig. 1-c) and the setting screw (Fig. 1-e). The ratchet head contains the ratchet wheel (Fig. 1-a) which is available in different variants depending on application.



Abb. 1

The ratchet head can be locked in order to lock the maximum torque (See Section 1.1). If the ratchet head is not locked a torque of between 10 – 30 Ncm can be set in an infinitely variable fashion via the setting screw (see Section 1.2).

Before initial use and following every use on the patient it must be disassembled, cleaned, dried and sterilised **before the next use** (see Section 3).

1.1 Setting the Ratchet Function

The ratchet head is supplied in a non-locked state. If the ratchet head is to be locked, turn the setting screw to less than 10 Ncm (Fig. 4) so that the spring is fully released.

Then pull the ratchet handle out of the ratchet head (Fig. 2), turn the ratchet handle through 90° in a clockwise direction and then engage the ratchet head (Fig. 3). After this set a torque of 20 Ncm so that the ratchet remains locked.

Attention: In the rigid state aware max. allowed load of 100 Ncm!

1.2 Setting the Torque Function

First ensure that the ratchet is not locked. Turn the setting screw until the required torque is achieved. The torque to be set is achieved when the corresponding indicator ring of the setting screw is still just visible on the end of the handle (Fig. 4). Under loading the ratchet head now clicks when the set torque is reached (Fig. 5). In no circumstances should it be turned further as otherwise there is a risk that the torque is significantly exceeded.

Important information:

If the ratchet is not used the setting screw is to be turned back so far until the spring is completely released.

1.3 Change of Tool Holder

(See 2.1-1 and 2.2-5)

2 Maintenance of the Torque Ratchet

2.1 Disassembly

1. Pulling both pin ends (red arrow) in the direction of the ratchet handle removes the ratchet wheel (Fig. 6-a) from the ratchet head (Fig. 6-b).
2. Unscrew the setting screw (Fig. 6-e) completely from the ratchet handle (Fig. 6-c).
3. Remove the ratchet handle (Fig. 6-c) from the ratchet head (Fig. 6-b) and remove the spring (Fig. 6-d).
4. Disinfection, cleaning and drying of the parts (see Section 3).
The ratchet is then assembled **before sterilisation**.


Abb. 6

2.2 Assembly

1. Before assembly all parts must be disinfected, cleaned and absolutely dry. Oil the points marked with a .
2. Insert the ratchet handle (Fig. 6-c) via the threaded rod onto the ratchet head (Fig. 6-b).
3. Insert the spring (Fig. 6-d) in the ratchet handle (Fig. 6-c).
4. Screw in the setting screw (Fig. 6-e) into the ratchet handle (Fig. 6-c) until the resistance of the spring is perceptible.
5. Pulling at both pin ends (red arrows) in the direction of the ratchet handle, sets the ratchet wheel (Fig. 6-a) into the ratchet head (Fig. 6-b). **Do not use force!**

3 General Preparation

3.1 Disinfection

Before initial use and immediately following every use the individual parts are to be disinfected in a suitable disinfection bath in accordance with the disinfection manufacturer's instructions.

Both chemical as well as thermal disinfection processes can be used which comply with the requirements of medical instrument disinfection.

3.2 Cleaning

Immediately following disinfection, the individual parts are to be thoroughly cleaned with a suitable cleaning agent in accordance with the cleaning agent manufacturer's instructions. Otherwise the pollutants dry on the parts and make cleaning more difficult.

In order to avoid corrosion and discolouring, the cleaning agent manufacturer's instructions regarding concentration, soaking time and temperature are to be strictly observed.

Coatings, pollutants etc. are to be thoroughly removed using a soft brush (never use metal brushes, abrasives, steel wool or excessive force).

Please note:

Observe the manufacturer's instructions for use when using an ultrasonic device for disinfection and cleaning.

3.3 Rinsing and Drying

Immediately after disinfection, cleaning or ultrasonic treatment, the parts must be carefully rinsed in clear, non-saline or distilled water (no mains water). In order to avoid water stains forming on the parts, they should be carefully dried afterwards. Dispensing with the rinsing process, not removing all residue and insufficient drying can lead to discoloration and traces of corrosion on the parts during the subsequent sterilisation.

3.4 Care and Inspection

Articulated parts must be treated with a sterilisable oil spray (paraffin-based) following the cleaning and disinfection stages. This spray reduces the friction of metal on metal, avoids fretting corrosion and reduces the risk of oxidation during heavy-duty use.

3.5 Sterilisation

Only carefully cleaned, dried parts displaying no signs of corrosion are permitted to be sterilised.

Poorly cleaned and insufficiently dried parts can lead to burning and corrosion damage.

Sterilising rusting parts with parts unaffected by rust can lead to the destruction of the other parts (flash rust).

Observe the instructions and recommendations of the manufacturer.

Please note:

The torque ratchet should be sterilised before using it for the first time.

Important information:

If the ratchet is not used the setting screw is to be turned back so far until the spring is completely released.